LACRIVERA

INSTRUCTIONS FOR USE

VeraPlug™

Sterile Preloaded Punctal Occlusion System

EN

Description

The VeraPlug™ punctal occluder is a medical device designed to occlude the punctum and thereby provide reduction or elimination of tear drainage through the inferior or superior puncta, thus maintaining lubricating tears on the surface of the eye. Each VeraPlug™ punctal occluder is molded from silicone. The VeraPlug™ is available in four sizes (small, medium, large and x-large). Each box includes two sterile occluders preloaded on an inserter, along with device identification labels for use by a qualified eye care physician, and instructions for use.

Indications for Use

The VeraPlug™ is for use in patients with dry eye syndromes.

Contraindications

Contraindications include, but are not limited to, eye infections, sensitivity or allergies to the occluder material and/or materials used in the manufacture of the device, blockage/infection of the lacrimal systems, inflammation of the eyelid, and epiphora.

Precautions

The VeraPlug™ may enhance the effect of some ocular medications in the eye. Depending on the type of medication being used the dose may need to be altered accordingly. Conditions such as blepharitis or ocular surface inflammation should be treated prior to use of punctal plugs. If the patient experiences irritation, infection or epiphora after the insertion of the VeraPlug™, the occluder should be removed.

Potential Adverse Events

The following complications may occur and the patient must be informed prior to use:

- Epiphora
- Pyogenic granuloma
- Foreign body sensation
- Infection of the lacrimal system
- Plug dislodgement or migration possibly requiring surgical intervention
- Punctal erosion
- Washout

Product Features & Documentation

Each box contains two individually packaged sterile VeraPlug™ punctal occluders preloaded on inserters for single use only. Additional accessories are not necessary for insertion. The use of stainless steel forceps may be recommended for removal. Instructions for use, patient implant card, and patient labels are also included in each box. A patient label should be applied to the implant card and provided to the patient.

Proper Sizing

The VeraPlug™ punctal occluder is available in four sizes: small, medium, large and x-large. To determine the proper size, begin with the smallest VeraPlug™ gauge (0.5mm) and insert the head of the gauge into the punctal opening. The gauge should fit snugly with a small amount of resistance. If the gauge enters the punctal opening with no resistance, the next larger size gauge should be tried in the same manner. Repeat this process until the proper size VeraPlug™ is determined based on the sizing chart shown.

Prior to Insertion

Prior to use, confirm the device is suitable for insertion. This includes visual confirmation that the plug is loaded properly on the inserter with no damage or debris and that, if applicable, the sterile barrier is intact by ensuring no damage or deformities are present that may compromise the sterile barrier of the encased device. Patients with intermittent tearing should receive probing and irrigation with a sterile irrigation solution to rule out pre-existing canalicular obstruction.

Insertion

- 1 Anesthetize the area of the punctum with a topical anesthetic placed in the conjunctival sac.
- 2 Apply a drop of saline solution or artificial tears onto the VeraPlug™ to help ease insertion.
- 3 Position the insertion instrument by placing the forefinger on the release button of the inserter and placing the occluder end of the insertion instrument over the patient's (superior or inferior) punctum.
- 4 Vertically insert the VeraPlug™ by positioning the occluder into the punctum until the cap is flush with the punctal opening. **FIGURE A**
- 5 When the occluder is properly seated, depress the release button and withdraw the insertion instrument. **FIGURE B**
- 6 Verify that the occluder is properly placed by confirming that the cap is flush with the punctal opening. **FIGURE C** After insertion, monitor the placement and integrity of the occluder to determine if/when the occluder may need to be replaced.

Removal

Should removal be indicated, grasp the vertical shaft of the occluder underneath the exposed cap with sterile forceps. Gently pull upward until the plug is removed. **FIGURE D**

Sterilization

VeraPlug™ sterile preloaded punctal occluders are offered in individual trays, two per box. The date of expiration should be confirmed prior to use. If the expiration date has lapsed the occluder should be discarded. Do not resterilize plugs as this may increase the possibility of post-operative infection.

Storage

The device should be stored in a dry place at 15-30° Celsius.

Warnings

The VeraPlug™ punctal occluder is intended for single use. Do not reuse. If the device is reused this may increase the risk of complications including, but not limited to, infection of the lacrimal system and plug dislodgement. If the sterile packaging is damaged or opened sterility is not guaranteed and the VeraPlug™ should be discarded. Any serious incident that has occurred should be reported to the manufacturer at info@lacri vera.com or by calling the manufacturer in the US at +1-855-857-0518. Serious events may also be reported to the competent authority in the country in which the user and/or patient is established. A Summary of Safety and Clinical Performance may be accessed by contacting the manufacturer. U.S. federal law restricts the sale of this device by or on the order of a physician.

Disposal


Any unused device and packaging may be disposed of through normal waste streams. In case of patient contact, the unused device should be disposed of as biomedical waste.

Device Materials

The device is composed of silicone with a titanium dioxide colorant. All materials are biocompatible. The device contains no materials of biological origin and no known endocrine disrupting devices nor carcinogens, mutagens, or reproductive toxins.

Authorized Representative

EC REP Emergo Europe, Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

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