

VeraPlug™ FlexFit™

Punctal Occlusion System

The VeraPlug™ FlexFit™ is designed for **simple sizing, easy insertion, patient comfort, and retention**—providing a proven treatment for chronic dry eye.

- Two sizes fit most patients for **simple sizing**
- FlexFit™ nose design provides **easy insertion** and **superior retention**
- Low-profile dome for excellent **patient comfort**

30-DAY RETENTION PLEDGE*

If a VeraPlug™ FlexFit™ is lost within 30 days of insertion, Lacriversa will replace it free of charge.



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
X-Small	VFF-5000
Small	VFF-5001
Medium	VFF-5002
Large	VFF-5003

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

SIZE	10 PAIRS
X-Small	BFF-5000
Small	BFF-5001
Medium	BFF-5002
Large	BFF-5003



VERAPLUG™ FLEXFIT™ SIZE CHART

PUNCTAL OPENING	SIZE
0.2mm to 0.3mm	X-Small
0.3mm to 0.5mm	Small
0.6mm to 0.8mm	Medium
0.9mm to 1.0mm	Large