

Dry Eye Diagnostics



BioGlo™

Fluorescein ophthalmic strips for evaluating contact lens fittings, corneal and conjunctival problems and applanation tonometry. CE marked. 100 or 300 individually wrapped sterile strips per box.



Lissamine Green

Ophthalmic test strips for evaluating corneal and conjunctival problems. Unlike Rose Bengal, it is sting-free. 100 individually wrapped sterile strips per box.



TearFlo™

Schirmer-type tear test for measuring and recording tear volume. CE marked. 100 individually wrapped sterile strips per box.

Instruments

Stephens Intracanalicular Plug Forceps S5-1059

For use in insertion of collagen and synthetic extended inserts.



Gland Expressor S5-1057

88mm forceps with closed paddles for use in the expression of the meibomian gland.



Stephens Meibomian Gland Compressor S5-1058

Overall length 105mm, paddle length 3.25mm. For use in the compression of the meibomian gland.



Sizing Gauge S8-1200 / S8-1205

0.5mm/0.6mm and 0.7mm/0.8mm for determining the proper size punctal occluder prior to insertion.



Miscellaneous

Dry Eye Compress LAC-1000

Featuring HydroHeat™ technology to provide clean, deep-penetrating, moist heat. Microwave for 30 seconds and apply 3-5 minutes a day to relieve dry eye and improve meibomian gland function, or freeze and apply cold to relieve swollen, puffy eyes.

