

VeraPlug™

Punctal Occlusion System

The **VeraPlug™** is designed to provide a simple and effective treatment for chronic dry eye.

- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**
- Available in four sizes



Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.

SIZE CHART

PUNCTAL OPENING	VERAPLUG™ SIZE
0.4mm to 0.6mm	Small
0.6mm to 0.7mm	Medium
0.7mm to 0.8mm	Large
0.8mm to 1.0mm	X-Large



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
Small	VSP-1001
Medium	VSP-1002
Large	VSP-1003
X-Large	VSP-1004

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

SIZE	10 PAIRS
Small	BNS-1001
Medium	BNS-1002
Large	BNS-1003
X-Large	BNS-1004